
CALIBREXTM organo 525

CALIBREXTM solutae 530

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SAFETY PRECAUTIONS

- Read operating instructions carefully, observe manufacturer's general recommendations and chemical compatibility limits.
- Refer to and follow regulations about handling of potentially hazardous reagents.
- Before each dosing operation, control working, tightness and bottle stability.
- Delivery jet should never point towards a person while dosing.
- Do not activate plunger with the screw-locking stopper attached and/or with flow control stopcock on STOP position.
- Seized parts should not be separated by applying force.
- Remove screw-locking stopper from delivery jet before dosing.
- When carrying dispenser or between usage, set volume on LOCK position, attach locking stopper and if available, set flow control stopcock on STOP position.
- Always hold bottle firmly when carrying the instrument.
- Do not remove delivery jet or flow control stopcock, unless instrument is properly purged.
- Prime and rinse dispenser carefully when left unused for a while or before disassembling.
- Take adequate measures to avoid electrostatic discharge when dosing flammable liquids.
- Operating temperature: 15°C - 40°C (60°F - 105°F).
- Instrument should only be used for its intended purpose.

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DESCRIPTION

Dispenser fits directly most laboratory bottles for safe, repeated dosing. Base thread 45 mm, supplied with three adapter of 32, 38 and 40 mm. Additional adapter sizes are available as options.

Built with materials carefully selected for their chemical compatibility; autoclavable fully assembled. CalibreX organo 525 includes a ground glass plunger, best suited for organics, non-crystallizing acid and base solutions.

CalibreX solutae 530 has a PFA coated plunger for weak or strong acids, bases and salt solutions. Each model available in volumes ranging 2.5-25 mL, 5-50 mL and 10-100 mL.

A flow control system with stopcock enables simple purge and recycling, thus limiting reagent waste. Calibration mechanism with integrated key located under plunger cap.

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DESCRIPTION

The dispenser is pre-arranged for being screwed directly on a flask for the distribution of reagents. A thread of 45 mm and 3 adapters of 32, 38 and 40 mm provided make the instrument compatible with the majority of flasks in the laboratory. Other sizes of adapters are available as an option.

Constructed with selected materials for their excellent chemical compatibility, each dispenser is autoclavable completely assembled.

CalibreX organo 525 possesses a glass rod piston for the distribution of organic solvents and non-saline solutions. CalibreX solutae 530 is provided with a piston with PFA coating which allows the dosing of acids, bases and other solutions which can crystallize. Each version covers the volume ranges of 2.5-25 ml, 5-50 ml and 10-100 ml.

A flow control system with stopcock makes purging and recycling of reagent simple and limits reagent waste. The calibration mechanism with integrated key is completely integrated.

OPERATION

Installation

Connect delivery jet to dispenser, screw by hand delivery jet adapter nut No. 15 (Fig. 1.1), without forcing. A longer delivery jet (150 mm) is available as an option; see accessories on page 4.

Cut feed tube No. 19 to desired length according to bottle size (Fig. 1.2) and insert on connecting body adapter (Fig. 1.3).

If needed, select an appropriate bottle adapter and install it on the instrument first. Use a stand to hold the unit if original reagent bottle does not provide for sufficient base stability; see accessories on page 4.

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Delivery jet position

Instrument rotates over 360° on connecting body (Fig. 2).

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DOSING INSTRUCTIONS

Volume setting


A dual graduation scale enables volume reading from a wide angle. Depress and move sliding cursor (No. 8) to desired volume (Fig. 3). Locking mechanism prevents unwanted move of sliding cursor to another position. A screw type button (No. 9) is supplied as an accessory with each dispenser. Follow Fig. 8 to change the sliding cursor. When using the screw type button, precisely align index to desired volume before locking, to avoid unwanted movement (Fig. 4).

Priming

Modell without stopcock

Remove stopper No. 17 and activate dispenser several times over short distance (≈ 30 mm) to eliminate any air bubble in the delivery jet. Collect liquid in a clean vessel for recycling (Fig. 5.1).

Modell with stopcock

Place selector in position  (Fig. 5.2) and remove stopper (No. 17). Slowly activate plunger. During distribution, the liquid will be directly recycled into the bottle. Repeat distribution several times over a short distance (≈ 30 mm) to eliminate any air bubble. Replace selector in distribution position . The instrument is then ready to use.

Dispensing

Remove stopper (No. 17) if not already done. Raising plunger head fills the dispenser barrel. Depress to dispense liquid. Best results are obtained by moving the plunger gently and regularly between upper and lower stops. Liquid density or viscosity may change dosing speed and instrument performance. Never distribute liquid if stopper is screwed on delivery jet or if the stopcock selector is in STOP position. Do not apply force when working.

End of dosing and storing

End dosing with plunger in lower position. Replace stopper (No. 17) on delivery jet and the stopcock selector on STOP position to prevent liquid leakage (Fig. 5.3). Purge and rinse instrument carefully before storing instrument.

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

Fig. 5.3

MAINTENANCE

Regular maintenance substantially contributes to instrument optimal working performance and long life span. Rinse when bottle has been emptied or before using another liquid.

Disassembling (Fig. 6)

Make sure dispenser has been purged before removing from the bottle. Pull feed tube No. 19.

- 1 Unscrew plunger head (No. 4) and re-move plunger. Never apply force to avoid any breakage.
- 2 Unscrew connecting body nut to access glass barrel, valve body (No. 13) and plate (No. 11). Do not misplace valve body ball (No. 12).
- 3 Unscrew delivery jet nut (No. 15) and separate from dispenser.

Cleaning

Soak dirty parts in a mild laboratory detergent. Take particular care of plunger, glass barrel and valve body. Residues on valve body seat or ball may lead to lack of tightness. A dirty valve body may block the instrument. Ultra-sonic bath may help remove residues.

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Reassembling

Damaged parts should be exchanged, using only original spare parts. To reassemble, proceed step by step in the disassembly reverse order. The exploded view on page 4 helps positioning each part. Make sure all connections are tight and check for correct aspiration and distribution functions of the instrument.

Calibration (Fig. 7)

Instrument is factory calibrated and delivered with an individual QC certificate.

Should performance of the instrument during routine control no longer correspond to manufacturer's specifications, or if conditions of use (liquid density, temperature, etc.) require recalibration, this can be easily performed using an analytical balance.

A new calibration is recommended when replacing either barrel/plunger assembly (No. 5), main body (No. 7) or valve body (No. 13).

- 1 Remove calibration safety seal (No. 1).
- 2 Remove cap (No. 2).
- 3 Lift calibration key (No. 3).
- 4 Adjust according to the +/- arrows engraved on calibration key. Each notch corresponds to a volume correction of:
525/530.025 ± 17.4 µL
525/530.050 ± 34.7 µL
525/530.100 ± 52.1 µL

Replace cap (No. 2) and apply a new calibration safety seal (No. 1).


Fig. 6

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
Fig. 8.1

Sterilization (Fig. 6)

Instrument has been designed for sterilization in the autoclave at 121°C / 250°F (20 min. at 1 bar), fully assembled. Do not sterilize an instrument which has not been purged and rinsed. Make sure delivery jet (Nos. 15, 16, 18) is correctly mounted.

Slightly untight connecting body nut (2), delivery jet nut or stopcock nut (3), as well as stopper (4). Place instrument horizontally in the autoclave, avoiding any direct contact with metal.

Instrument must be dry and cooled down to room temperature (2 hours) before securely re-tightening all connections, without excessive force.

It is recommended to verify inaccuracy after the first few autoclaving cycles in order to guarantee instrument performance.

Repeated autoclaving may affect material colours. Correct autoclaving and resulting sterility are the responsibility of the user.

Materials

All parts in contact with the liquid are chemically inert:

- Feed tube: FEP
- Valve body: Ceramic
- Valve ball: Ceramic
- Valve spring: Platinium-iridium
- Plate: PTFE
- Barrel: Borosilicate glass
- Plunger: Ground glass (525) or PFA coated glass (530)
- Connecting body: ETFE
- Delivery jet: FEP/PCTFC
- Stopper: ETFE

Warning: Vapors of concentrated strong acids may damage your instrument after some time and reduce its life span. It is recommended not to leave the dispenser permanently on a bottle containing concentrated strong acids and to rinse it after use.

To prevent long term loss of product by evaporation, the same recommendation is valid when using volatile solvents. In case of doubt, check for chemical compatibility or refer to manufacturer before use.

The QR code printed on the instrument provides access to the latest chemical compatibility chart.

Damages due to incompatible reagents are not covered by the warranty.

USE OF ACCESSORIES

Screw type button (Fig. 8)

To change from the sliding cursor to screw type button:

- 1 Unscrew plunger head (No. 4)
- 2 Remove plunger.
- 3 Make sure body is in the upper position. Remove upper apron (No. 6) by lifting front side.
- 4 Slide cursor toward the top and remove it from the body.
- 5 Insert screw type button, slide it toward the bottom and screw it to lock.

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Replace upper apron (No. 6) until it "clicks" (Fig. 8.1)
Insert plunger and securely tighten. A calibration is not necessary after this operation.
Note: Upon reinstalling sliding cursor, make sure notches are pointing down.

Long delivery jet

If required, a long delivery jet (150 mm) can be ordered as an accessory. Screw stopper (No. 17) on delivery jet. Unscrew delivery jet adapter from the dispenser to install replacement.
The addition of the stopcock extends overall delivery jet length by 20 mm.

Sterile air filter

Air inlet hole in the back of connecting body can be fitted with a membrane filter with Luer connection. Use a sharp blade to cut out air inlet protection cover and insert filter in the opening. Filter must be removed before autoclaving the instrument.

WARRANTY

Dispensers are guaranteed for a period of two years against any factory or material defect. Glass breakage or cases due to non respect of manufacturer's and safety instructions are excluded from the warranty. Repairs and changing parts do not extend the warranty period.
Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased, after obtaining return authorization.
Purge and decontaminate the instrument prior to returning it .

PERFORMANCE DATA AND ORDERING INFORMATION


| Volume mL | Increment mL | Inaccuracy (E%) | | | Imprecision (CV%) | | | Calibrex 525 Cat. No. | Calibrex 530 Cat. No. |
|-----------|--------------|-----------------|-----------|-----------|-------------------|-----------|-----------|-----------------------|-----------------------|
| | | Min. vol. | Mid. vol. | Max. vol. | Min. vol. | Mid. vol. | Max. vol. | | |
| 2.5 - 25 | 0.5 | <±/- 1.5% | <±/- 1.1% | <±/- 0.6% | < 0.5% | < 0.35% | < 0.1% | 525.025 | 530.025 |
| 5 - 50 | 1.0 | <±/- 1.5% | <±/- 1.1% | <±/- 0.6% | < 0.5% | < 0.35% | < 0.1% | 525.050 | 530.050 |
| 10 - 100 | 1.0 | <±/- 1.5% | <±/- 1.1% | <±/- 0.8% | < 0.5% | < 0.35% | < 0.1% | 525.100 | 530.100 |

Performance values obtained by a regular smooth and steady movement, with bidest, water at constant temperature (±0.5°C) comprised between 20 and 25°C, in accordance with ISO 8655.
Omission to tighten connecting body ring before autoclaving, and/or over-tightening ring when dispensing, may reduce performance.

ACCESSORIES

| Cat. No. | Description | Packaging |
|-----------|--|-----------|
| 1.525.FC | Stopcock, all models | 1/pack |
| 1.525.120 | Standard delivery jet, 25mL | 1/pack |
| 1.525.121 | Standard delivery jet, 50mL/100mL | 1/pack |
| 1.525.150 | Long delivery jet, 25mL (150 mm) | 1/pack |
| 1.525.151 | Long delivery jet, 50mL/100mL (150 mm) | 1/pack |
| 525.350 | Replacement feed tube, 350 mm | 1/pack |
| 525.706 | Feed tube by the meter, uncut | meter |
| 320.SB100 | Dispenser stand, fits all models | 1/pack |
| 320.BC100 | Stand for remote liquid intake* | 1/pack |
| 1.525.525 | Calibration safety seal | 25/pack |
| 1.525.526 | Calibration safety seal | 100/pack |

*Dispenser and feed tube to be ordered separately.

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DONNEES TECHNIQUES ET CATALOGUE

| Volume mL | Max. cinematic viscosity (cst) | Example of liquid | Test conditions for all models (without stopcock) |
|-----------|--------------------------------|-------------------|---|
| 25 | 400 | Silicone oil | Dispensing time <10s. |
| 50 | 250 | Silicone oil | Dispensing force <3kg. |
| 100 | 150 | Silicone oil | Pressure in barrel <3bar |

Cinematic viscosity is indicated in centistokes (cst). The dynamic viscosity in centipoises (cps) is obtained by multiplying the cinematic viscosity by the density of the liquid (kg/m³).

ZUBEHÖR

| Cat. No. | Description | Packaging |
|------------|--|-----------|
| 1.525.GLxx | Threaded adapter (PP) ∅ 32, 38, 40 mm | 1/pack |
| GLP.xx | Nested threaded adapter (PP) to be used with 1.525.GLxx ∅ 22, 25, 28, 30, 34, 36 | 1/pack |
| 314.xxxx | Amber glass reservoir (∅ 32/45 mm) 0250, 0500, 1000, 2500 mL | 1/pack |
| 314.xxxxPE | Amber glass reservoir (∅ 32/45 mm) 0500, 1000, 2500 mL with anti-spill PE coating | 1/pack |
| 314.0500P | Pyrex glass reservoir (∅ 32 mm) 500 mL with side neck for filter | 1/pack |
| 315.xxxx | Polyethylene reservoir (∅ 32/45 mm) 0250, 0500, 1000, 2500 mL | 1/pack |

Replace xx by adapter diameter and xxxx by reservoir volume.

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Chemical resistance
Scan QR code to access chemical resistance chart.