# SOCOREX

## MICROPIPETTE 841

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Your new solid interface pipette has two essential features: it is simple and accurate. By following the few recommendations of this leaflet, you will take full advantage of your instrument.

#### DESCRIPTION

This is an adjustable pipette with five pre-set volumes. The plunger draws the sample into the precision glass capillary without any dead volume. No air interface means that the volume dispensed will always be the same regardless of viscosity (ideal for serum). The plunger is designed to deliver the sample without any noticeable residue, which allows the capillary to be re-used with other fluids.

5 pipette sizes cover the full range of usual sample volumes from 1 to 200 ul.

Each pipette is individually computer-tested and supplied with a print-out.

Volume selection is a matter of seconds. Simple mechanical device gives high precision, cuts out risk of error. At all times, the selected volume clearly appears on the barrel.

#### INSTALLING GLASS CAPILLARY

Set pipette on **largest volume**; loosen collet A. Slip capillary over plunger **with colour mark first** until it bottoms against stop.

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


## CALIBRATION

Line up calibration mark of capillary with flat end of plunger, as shown in enlargement fig. 1. Tighten collet. The pipette is ready. When plunger rod is completely depressed; tip will protrude slightly from glass tube.

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## VOLUME SELECTION

Barrel C will slide after rotating it from locked position ▼ to selector position I. Place selected graduation mark, for example 50 ul, flush with end of ring D; lock into position ▼ (fig. 2).

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## REPLACING GLASS CAPILLARY

Loosen collet A and remove used capillary. Set pipette on **largest volume** and insert new capillary, then proceed as above. **Do not try to remove a broken capillary over the plunger;** this would damage it. Remove plunger rod and remove broken glass over the rod. Make sure that no glass has been left in the collet A.

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## COLOUR CODE


The colour appearing on the finger-rest of the pipette is to be found on the matching capillaries and spare plungers labelling. Only assemble components of the same colour.

## PIPETTING

1. Depress plunger fully and immerse capillary into solution. Release slowly. Plunger will glide freely once it is wet.
2. Wipe capillary being careful not to touch opening.
3. Place capillary against wall of receiving vessel and smoothly depress plunger until it butts against stop.
4. After a series of tests, discard capillary or rinse in distilled water or saline by operating plunger rapidly a few times.

## REPLACING PLUNGER

To remove a worn plunger, lift cap E off finger button and use built-in Allen wrench to loosen set screw F. Turn pipette upside and insert a new plunger until it bottoms. Tighten set screw F and press cap E on button.

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