
# CALIBRA<sup>®</sup> 852

## DIGITAL MULTICHANNEL PIPETTE DIGITAL MEHRKANALPIPETTE PIPETTE DIGITALE MULTICANAUX

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#### DESCRIPTION

Adjustable volume multichannel micropipette for use with detachable polypropylene tips, for an easy and continuous work with standard microplates.

The sample is drawn into the polypropylene tip without any contact with the plunger assembly. The integrated tip ejector is located on the handpiece allowing the automatic release of the tips.

Simultaneous delivery of 8 or 12 samples from 1 µl to 200 µl by means of three pipette sizes, autoclavable at 121°C fully assembled.

Each instrument is individually quality controlled and comes along with its calibration certificate.

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#### HOLDING THE MULTICHANNEL MICROPIPETTE (Fig. 2)

The ergonomic shape of the Calibra<sup>®</sup> greatly reduces operator fatigue. Place fingerrest on the phalanx of forefinger with ejector knob turned towards yourself. Ease of access of both push button and ejector knob allows natural and comfortable thumb movement.

The multichannel casing revolves on the handpiece and allows the selection of the best working position.

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#### VOLUME SETTING (Fig. 1)

Two cylindrical cams fitted with pre-calibrated steps rotate on their axis. The large cam shows the figures on the left hand side of the window, the small one on its right hand side.

Setting the volume as described below and illustrated in Fig. 1 is fast, easy and glove friendly.

1. Push plunger button with thumb and pointer to the first stop and lock it by a clockwise quarter of a turn.
- 2a. Turn setting wheel clockwise or counter-clockwise until the expected key volumes appear at the left of the window.
- 2b. If required, pull out setting wheel (position clicks-in) and repeat point 2a for the fine tuning.
3. Push setting wheel back to the initial position and unlock plunger button.

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Abb. 1. Volumen einstellen

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## TIPS AND FILTERS

Disposable tips with or without filter were carefully designed to provide the best pipette and tip combination for all types of applications. See p. 4 for ordering information. Tips are supplied on 96 piece racks, with the 9 mm spacing fitting the instrument channels.

Place the instrument on a row of tips, press firmly to ensure a solid fit of the tips on each channel to avoid leakage.

Rinsing every new tip once or more with the pipetting solution (fill – empty – expulse) will improve precision.

Regular tips are autoclavable at 121°C. Filter tips are supplied in sterile and autoclavable versions.

## TIGHTNESS

To work correctly, the air chamber of each barrel has to be tight. Leaks can be detected by observing the tip: a droplet forming at the extremity signals a leak in the channel concerned.

Following spots have to be watched for leakage:

1. Matching of the tips on nozzles
2. O-rings providing tightness to the plungers
3. Damaged nozzles

Leakage will usually be eliminated by changing the tip. If the problem remains, change the suspicious O-ring. The wear of the O-ring is very slow under normal use, somewhat faster under frequent autoclaving. Sterilization therefore requires higher frequency tightness control. Damaged nozzles should be changed.


Fig. 2

## PIPETTING OPERATION (Fig. 3)

Forward pipetting mode:

**Phase 1** Depress plunger button to the intermediate stop ending calibrated stroke before immersing tips into appropriate reagent container.

**Phase 2** Immerse 8 or 12 tips vertically 5 mm into liquid simultaneously and smoothly release plunger button. Withdraw tips from liquid vertically after 2-3 seconds, sliding them along reservoir wall.

**Phase 3** Place tips against wall of microtiterplate, slowly depress plunger button to first stop and let liquid out. Depress button to second stop: the excess air volume produced by depressing the plunger beyond the calibrated stroke blows out possible fluid particles which may cling to the tips. Remove the multichannel pipette from the wells of the microplate, with plunger knob in its lowest position. Release knob slowly.

**Tip ejection** To eject used tips, press ejector knob E, as shown in above illustration.

Apply new tips for every liquid; change tips if droplets are retained after blow out.

Reverse pipetting mode also possible with the Calibra® multichannel pipette.

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## CALIBRATION

Each Calibra® has been calibrated permanently at the factory and individually controlled. The Q.C. certificate accompanying the instrument includes control data, serial number and operator identification. Calibration is carried out at a constant temperature (+/- 0.5°C) comprised between 20 and 25°C with Aqua bidest, according to ISO 8655.

Under normal conditions, a Solid Calibration® pipette does not need to be recalibrated. After replacing a part, which could influence the measure, the calibration should be checked, preferably by gravimetry according to the above mentioned norm. Ten measurements should be performed after checking the seal tightness of each channel of the instrument.

**In-lab calibration (Fig. 4):** Whenever required, for instance when pipetting viscous or non-aqueous liquids, calibration adjustments can be performed by laboratory personnel. Reach calibration screw c) located inside plunger button by removing comfort pad a). By means of M3 allen key b), delicately turn calibration screw clockwise to decrease volume, and vice versa, using values indicated below. Check resulting volume after each fine tuning. Cover push-button with comfort pad again.

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## MAINTENANCE

The Calibra® 852 pipettes are in conformity with Council directive IVD 98/79 EEC (in vitro diagnostic medical devices). Consequently, any change on the pipettes (i.e. parts, calibration) is done under responsibility of the user. The models were designed to sustain long-term operation with minimal maintenance and service requirements.

After a longer period of work, it may become necessary to change an O-ring (droplet forming at the tip) or to clean the plungers (plungers dragging).

Changing the O-rings: to have access to the O-rings No. 2 and 9 (lower assembly), remove the two lock pins No. 62 and pull barrel assembly out of casing. Free ejector head No. 3 from claws of ejector shafts No. 12, remove circlips No. 63 and pull plungers out of barrels. Remove barrel holder No. 10 held in place by upper circlip No. 65. Barrels and O-rings are now free for inspection.

Slightly grease O-rings before re-assembling unit.

Cleaning the plungers: clean the plungers with a mild detergent, rinse, add some grease and wipe with a clean cloth.

Clean, grease and check instrument once a year or more depending upon frequency of use.

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## STERILIZATION

The unit has been designed for repeated sterilization in the autoclave at 121°C (20 minutes), fully assembled. Before autoclaving, set pipette to the largest possible volume. Place instrument horizontally in the autoclave avoiding any direct contact with metal.

Allow instrument to dry and cool down before use. Tighten lower assembly nut if loose. If plunger movement is not smooth, open the pipette and slightly lubricate plunger with silicone grease supplied with the instrument.

Check tightness and performance after first autoclaving cycle then regularly but at least after 50 autoclaving cycles.

A slight change of material colours may appear after long term repeated autoclaving.

Correct autoclaving and resulting sterility are the responsibility of the user.

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## WARRANTY

Your Calibra® 852 pipette is guaranteed against any material or manufacturing defects for the period of time specified in its QC certificate. Damages due to non-respect of manufacturer's instructions, safety precautions or autoclaving conditions, as well as material colour alteration are excluded from the warranty. Repair and replacement of parts do not extend warranty time. Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased after obtaining return authorisation. **Decontaminate instrument prior to returning it.**

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**PERFORMANCE DATA AND ORDERING INFORMATION**

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

A = Autoclavable, R = Rack, B = Bag, L = Long, F = Filter, S = Single wrapped

Accuracy* (mean error %)	Precision* (CV %)
1 µl < ± 3,5%	1 µl < 3,0%
10 µl < ± 1,5%	10 µl < 1,0%
10 µl < ± 1,0%	10 µl < 1,0%
100 µl < ± 0,7%	100 µl < 0,4%
20 µl < ± 0,9%	20 µl < 0,6%
200 µl < ± 0,7%	200 µl < 0,3%

\* Above values are actual Socorex Q.C. limits achievable in daily routine with high quality tips. (Aqua bidest., temp. 20-25°C (± 0.5°C) per ISO 8655)

Description	Capacity	Size	Cat. No.	Description	Capacity	Size	Cat. No.	Description	Cat. No.
Single well reservoir	70 ml	20/Pack	<b>330.01</b>	8 channel reservoir sterile	40 ml	10/Pack, single wrapped	<b>330.08.9</b>	Work station holding 3 multichannel pipettes	<b>320.340</b>
Large reservoir with lid	175 ml	20/Pack	<b>330.01 L</b>	12 channel reservoir sterile	40 ml	10/Pack, single wrapped	<b>330.12.9</b>		

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