
SOCOREX® *ultra* 1810

ANIMAL HEALTH SYRINGE VETERINÄR SPRITZEN SERINGUES VÉTÉRINAIRES

OPERATING INSTRUCTIONS BEDIENUNGSANLEITUNG MODE D'EMPLOI

SAFETY PRECAUTIONS

- Read the operating instructions carefully before use. Keep available for future reference
- Observe the manufacturer's recommendations
- Follow instructions for proper use of injectables
- Before each use, check instrument for perfect working conditions
- Purge, rinse, dry syringe and lubricate the tightness parts carefully after each use
- Never try to disassemble blocked parts and do not use force to separate plunger from cylinder or valves
- Repeat use of the same needle may carry over contaminations to other animals
- Instrument manufacturer cannot be made responsible for needle breakage
- Syringe should be used by professionals only
- Mind possible hand fatigue and its medical consequences during prolonged vaccine campaign
- Do not use syringe for purposes other than intended for
- A saline environment may hasten replacement of instruments parts
- Do not force on handle if needle is blocked. High pressure inside syringe, may lead to glass cylinder breakage and valve damages
- Only use original spare parts

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DESCRIPTION

Congratulations for buying a high quality Swiss product.

The automatic syringe *ultra* 1810 allows trouble free injection of volumes as small as 0.02 mL of various water or oil based vaccines with excellent reproductibility. The syringe with spring loaded plunger and three-way valve system is intended for micro-range injections (up to 0.1, 0.2, 0.3 and 0.5 mL) such as for fishes, chicks and small laboratory animals. Larger volumes (1 and 2 mL) ideal for poultry or pigs.

It can be autoclaved fully assembled at 121°C (20 minutes). Untighten valves and casing nut prior autoclaving.

Valves are self-contained, avoiding lost of small parts. Volume is easily set with a micrometer screw and locked by nut. The graduated glass barrel is interchangeable. Pistol handle syringe is supplied with feed tubing.

Veterinary needles and drenchers to be ordered separately. Ask for more details about our large selection of Luer Lock stainless steel needles.

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HOLDING THE SYRINGE

Ergonomic, user friendly handle, perfect balance and smooth plunger movement, ensure fatigue-free operation even during mass vaccination at rates of up to 1000 injections per hour.

Place hand on handle and move plunger with hand palm.

To change position of syringe casing (12), untighten screw (7), pull toward the body, rotate the syringe body into one of the eight possible slots. Push back body and retighten screw (7) keeping plunger in forward position.

OPERATION

The instrument is supplied ready to use, no assembling required.

- Aspiration from bottle at belt: puncture bottle cap with feed needle and vent cannula and connect feed tubing to the feed needle. Filling by gravity the entire silicone tubing, facilitate priming. Verify absence of rubber residue in tubing before connection. Repeated puncture of the rubber seal is not recommended.
- Aspiration from bottle hanging from holder (neck toward the bottom): connect feed tubing to the optional draw-off needle with air inlet and filter (ref. 1.187.590).
- Aspiration from open bottle: connect sinker to the end of the feed tubing and immerse into liquid.

Priming

Set instrument on maximum volume. Prime feed tubing by activating plunger, until liquid completely fills the syringe. Hold syringe upright to expel all air bubbles and/or residues by activating plunger to minimum of five times. The lost liquid during priming can be collected in a clean recipient and returned to container. After changing bottle, verify proper priming (no air bubble) to guarantee optimal performances.

Volume setting (Fig 1 and 2)

Activate plunger a few times before setting volume to enable free plunger movement (stick slip). Loosen locking nut (5) and turn socket (4) to center O-ring on desired graduation on models 0.1, 0.2 and 0.3 mL (Fig 1). For 0.5, 1 and 2 mL, center plunger edge on desired graduation (Fig 2). Tighten locking nut (5) in order to prevent any volume alteration. Activate a few times the plunger to verify good operation.

Injection

Attach injection needle (to be ordered separately) to nozzle (mind adequate needle size). During injection, do not apply force on plunger rod and check output liquid from needle.

Too much pressure inside barrel can cause barrel breakage or leakage (risk of injury).


Fig. 1

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
Fig. 2

End of dosing

When dosing is completed, remove feed tubing and activate plunger several times to eliminate any residual liquid. Rinse, disassemble, clean and dry the instrument parts according to the instructions in the "maintenance" chapter.

Never store a syringe without prior purging, cleaning and drying parts. Lubrification of tightness parts (plunger O-ring, barrel and valves) avoids blockage over time. Medical grade white and silicone oils recommended. Other stable food oils can also be used.

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MAINTENANCE

Cleaning

Syringe and accessories must be carefully cleaned after each use.

Aspirate water or an appropriate cleaning solution (detergent) and rinse thoroughly. Solution should not contain any residue. Purge cleaning solution by activating plunger several time..

Disassembling

Remove needle and feed tubing. Hold handle (6) and unscrew casing nut (7) while holding syringe body (fingers on barrel).

Slowly pull forward casing to disengage plunger (8). Never apply excessive force if the plunger is blocked. Remove barrel, watch out not losing barrel washer (11). Untighten both connectors (15 and 16), remove valves (13 and 14). Rinse all parts with appropriate detergent or warm water (no residue). It is mandatory to let dry all parts (valves, connectors, syringe casing, plunger and barrel) before reassembling.

Valves reassembling (Fig 3)

Introduce each valve in its specific housing, respecting mechanical coding. Lubricate valve if syringe is stored. Screw inlet and outlet valves. Do not use tool, neither apply excessive force while manually tightening.

Plunger O-ring replacement (Fig 4 and 5)

Apply only to 0.1, 0.2 and 0.3 mL models. Use thin needle to remove O-ring. Avoid damaging the O-ring groove.

Mount new lubricated O-ring on guide (1). Push O-ring into the plunger groove (2).

Plunger reassembly (Fig 6 and 7)

Slightly lubricate barrel entrance (1), insert half of plunger into barrel without force (2). Place barrel into the syringe casing. Couple plunger onto plunger rod (3) and push the plunger fully inside barrel (4).

Introduce the syringe casing into the desired handle slot and partially screw casing nut (7). Position the barrel as to see the graduation. Push the plunger fully into the barrel by pressing on the piston rod (1) and tighten the nut (7) without excessive force. Assure correct syringe operation with several plunger movements. Plunger should move freely and clearly stop at end of travel.


Fig. 3

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
Fig. 6


Fig. 7

TROUBLE	POSSIBLE CAUSE	REMEDY
Erratic operation	Syringue used without liquid inside Plunger, O-ring and barrel not properly lubricated	Use with liquid Use appropriate lubricant
Blocked valve	Syringe not cleaned Valve parts blocked by residue	Disassemble and clean Soak valve and clean (remove residue)
Leaking valve	Residue in valve (i.e. rubber residue from septum, cristalization) prevent tightness Damaged valve	Clean valve Replace part
Wrong volume	See leaking valve	See leaking valve
Inadequate performance	Priming not completed, air bubble in system Barrel washer not positioned correctly or missing	Prime properly Change barrel washer
Liquid accumulation behind plunger	0.1 mL – 0.3 mL models, damaged O-ring or missing 0.5, 1 and 2 mL models, this may occur after long injection series	Check plunger, change O-ring Remove liquid when doing maintenance
Plunger blocked	Syringe stored without cleaning Needle clogged	Do not apply force to unblock. Disassemble and soak syringe in warm water and separate plunger parts Change needle
No priming possible	Check for blocked valve Check for leaking valve Check for clogged needle	Clean and reassemble Change valve if required Change needle if required

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STERILIZATION

The *Ultra* 1810 syringe can be autoclaved fully assembled at 121°C, during 20 minutes in a sterilization pouch (Qualitix cat. No. 3210.1326) laying directly on a metal grid. Do not sterilize an instrument that has not been previously purged and cleaned. See chapter "Maintenance".

Lubricate plunger O-ring, cylinder and valves before autoclaving the syringe to avoid blockage. See chapter "End of dosing".

Slightly untighten both valve connectors (15 and 16) and casing nut (7).

After sterilization, the instrument must be cool and completely dry prior to retighten all connectors and nuts.

Correct autoclaving procedure is the sole user responsibility.

MATERIALS


All parts in contact with liquid are made of carefully selected and chemically inert materials:

- Feed tubing: silicone
- Barrel: sodocalcic AR glass (0.5, 1, 2 mL) and borosilicate glass 3.3 (0.1 to 0.3 mL)
- Plunger: stainless steel
- Barrel washer, plunger O-ring: FPM (FDA compliant)
- Valve: PP, stainless steel and FPM (FDA compliant)
- Syringe casing: PA
 - Sinkers and feed cannula: nickel/chromed plated brass

PERFORMANCE AND ORDERING INFORMATION

Package contents

- Instrument
- Feed tubing (17)
- Sinkers (18)
- Feed needle (19) and vent cannula (20)
- Operating instructions
- Spare parts bag with barrel O-ring, plunger O-ring (except for 0.5, 1 and 2 mL models), guide for O-ring

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ZUBEHÖR

Cat. No.	Description	Packaging
Valves		
1810.915	Intake valve (0.1 - 0.5 mL)	1 pc
1810.945	Outlet valve (0.1 - 0.5 mL)	1 pc
1810.916	Intake valve (1 - 2 mL)	1 pc
1810.946	Outlet valve (1 - 2 mL)	1 pc
Spare part kit		
1.1810.72001	0.1 mL	1 /pk
1.1810.72002	0.2 mL	1 /pk
1.1810.72003	0.3 mL	1 /pk
1.1810.72005	0.5 mL	1 /pk
1.1810.7201	1 mL	1 /pk
1.1810.7202	2 mL	1 /pk
Various syringe accessories		
1.187.590	Draw off needle	6 /pk
Drenching cannulas		
376.0925S	Straight, 25 mm	3 /pk
376.1251S	Straight, 51 mm	3 /pk
376.16102S	Straight, 102 mm	2 /pk
376.0925B	Bent, 25 mm	3 /pk
376.1251B	Bent, 51 mm	3 /pk
376.16102B	Bent, 102 mm	2 /pk

ACCESSOIRES

Cat. No.	Description	Packaging
Premium lubricant ABF		
1.170.705	70 mL	1 pc
1.170.906	70 mL (pack of 6)	6 /pk
1.170.755	500 mL	1 pc
1.170.956	500 mL (pack of 6)	6 /pk
1.1810.700	White oil	1 /pc
Superior quality needles, Luer Lock *		
371.0503	0.5 x 3 mm / 25G x 1/8"	12/pk
371.0504	0.5 x 4 mm / 25G x 1/6"	12/pk
371.0603	0.6 x 3 mm / 23G x 1/8"	12/pk
371.0605	0.6 x 5 mm / 23G x 3/16"	12/pk
371.0608	0.6 x 8 mm / 23G x 1/3"	12/pk
371.0704	0.7 x 4 mm / 22G x 1/6"	12/pk
371.0706	0.7 x 6 mm / 22G x 1/4"	12/pk
371.0710	0.7 x 10 mm / 22G x 3/8"	12/pk
371.0806	0.8 x 6 mm / 21G x 1/4"	12/pk
371.0810	0.8 x 10 mm / 21G x 3/8"	12/pk
371.1010	1.0 x 10 mm / 19G x 3/8"	12/pk
371.1213	1.2 x 13 mm / 18G x 1/2"	12/pk
371.1420	1.4 x 20 mm / 17G x 3/4"	12/pk
371.1625	1.6 x 25 mm / 16G x 1"	12/pk

* Additional sizes, see www.socorex.com

Viscosity data

Volume mL	Maximal cinematic viscosity * (cst)	Example of liquids
0.1 to 0.5	200	Glycerol 90%
1	1200	Silicon oil
2	400	Glycerol 94%

Test conditions to determine maximum viscosity

Syringe with horizontal 1 meter feed tubing. Syringe set at maximal volume. Aspiration time < 3 sec. Needle diameter 0.8mm (21G).

* Cinematic viscosity in centistokes (cst). Dynamic viscosity in centipoises (cps) obtained by multiplying cinematic viscosity by liquid density (kg/dm³). Performance data and viscosity limit are related to selected needle size.

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Volumen mL	Maximal cinematic viscosity * (cst)	Beispiel an Flüssigkeiten
0.1 to 0.5	200	Glycerol 90%
1	1200	Silicon oil
2	400	Glycerol 94%

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WARRANTY

The *ultra* 1810 syringe is guaranteed for a period of one year against any manufacturing or material defect. Cases due to non-respect of manufacturer's instructions are excluded from the warranty. Glass barrel breakage, damages to the plunger or valves due to improper handling or cleaning are not covered by the warranty.

Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased, after obtaining a return authorisation.

Rinse, dry and decontaminate instrument prior to its return.

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