

Acura[®] *electro*



your electronic choice









Operating instructions Models 925, 935 and 955

4.8V version



 **SOCOREX**
SWISS











TABLE OF CONTENTS

| | |
|---|--------------|
| Foreword | 3 |
| Safety precautions and application limits | 4 |
| Pipette handling | |
| Maintenance and sterilisation | |
| Charging stand, power handle and power supply | |
| Instrument description | 5 |
| View | |
| Keys and functions | |
| Start button | |
| LCD display | |
| Power handle | |
| Optional charging stands | |
| Power supply | |
| Use of the instrument | 6 |
| Supply contents | |
| Inserting power handle | |
| Replacing power handle | |
| Charging power handle | |
| Programming | 7-16 |
| Holding the electronic pipette | 7 |
| Left or right display reading | 7 |
| Forward mode, programming and pipetting | 8 |
| Reverse mode, programming and pipetting | 9 |
| Stepper mode, programming and pipetting | 10-11 |
| Stepper mode, excess volume control | 11 |
| Stepper mode, interruption of pipetting sequence | 12 |
| Dilution mode, programming and pipetting | 13-14 |
| Dilution mode, interruption of pipetting sequence | 15 |
| Mixing | 16 |
| Operation | 17-21 |
| Selection of pipetting speed | 17 |
| Charge level of power handle | 17 |
| Changing volumetric module | 18-19 |
| Programming for volumetric module | 20 |
| Correction of volumetric module selection | 20-21 |
| Use of Pasteur pipette | 21 |
| Adjustment of tip ejector | 21 |
| Maintenance and sterilisation | 22-23 |
| Cleaning | 22 |
| Replacement of tightness parts | 22-23 |
| Sterilisation | 23 |
| Instrument calibration | 24-25 |
| Calibration increments | 24 |
| New calibration | 24-25 |
| Trouble shooting | 26-27 |
| Error messages | 26 |
| Instrument reset | 26-27 |
| Other failures | 27 |
| Performance | 28 |
| Warranty | 29 |
| Ordering information | 30-31 |
| Instruments, accessories and tips | |
| Exploded drawings | IV-VI |

Foreword

Congratulation for choosing a Socorex instrument. You purchased a superior quality product built to last, and to satisfy you for a long time.

The Acura® *electro* pipette allows precision liquid handling without hand fatigue. The microprocessor-controlled instrument includes an extended life NiMH battery pack for high performance.

Main user benefits

- Ergonomic, optimized weight and hand balance for perfect working comfort
- User-friendly through intuitive software programme, easy to use
- Display selection for left or right-handed user
- Extended pipetting autonomy, instant exchange of battery pack eliminates interruption of work
- Autoclavable volumetric modules are interchangeable on single control unit
- Instrument software covers all volume modules
- Justip™ ejector facilitates fitting and ejection of tips

Before using the instrument for the first time, read instructions carefully. Special attention should be paid to the safety precautions and application limits. Keep this booklet for future reference.

SAFETY PRECAUTIONS AND APPLICATION LIMITS

Before using the instrument for the first time, read the safety precautions and application limits carefully. Socorex will not assume any responsibility for problems related to erroneous use of the instrument.

Pipette handling

- **Use only 4.8V power handle (battery) with appropriate 4.8V devices**
- Refer to and follow regulations about handling of potentially hazardous reagents
- Before use, check tip/nozzle tightness and working condition of the instrument
- Emergency stop: press the “set/stop” key to interrupt immediately any plunger movement in stepper and dilution modes
- Do not use the Acura[®] *electro* in areas where there is explosion or flame hazard
- Do not place the Acura[®] *electro* on charging stand with a filled tip
- Do not let liquid penetrate inside the pipette housing (upper assembly)
- Change volumetric module **only** with charged power handle in place
- Do not use any volumetric module (lower assembly) not specifically intended for your instrument
- Instrument should not be used at temperatures below 5°C and over 40°C
- This product should be used only for its intended purpose
- Mind possible hand-fatigue during serial pipetting and its potential medical consequences (such as repetitive strain injuries, RSI)

Maintenance and sterilisation

- Do not use aggressive solutions (such as acetone) to clean volumetric modules. Use water or alcohol (ethanol) instead
- Only the volumetric module (lower assembly) is autoclavable at 121°C. The control unit (upper assembly) is **not** autoclavable
- No liquid must penetrate inside the pipette control unit (upper assembly) nor the changing stand
- Refer to instructions for any operation (maintenance, change of volumetric modules) carried out on the instrument
- Servicing of control unit must be done by authorised, trained personnel only
- Use only Socorex original parts and accessories (power handle, pipette tips, spare parts)

Charging stand, power handle and power supply

- **Use only 4.8V power handle (battery) with appropriate 4.8V charging stand**
- Power handle (battery) must be completely charged before first use
- Power handle can be charged either separately, or when placed on the instrument
- For best power handle durability, avoid recharging when "low bat" is not lit
- Do not use any power supply other than the original one supplied by manufacturer
- Do not damage electrical cord or plug of charger with heavy or sharp items
- Do not expose pipette, power handle, charging stand or power supply to heat or liquid spillage
- Remove power handle from the instrument if it is stored for a long period of time; this prevents discharge of the batteries
- Recycle used or damaged power handle according to applicable local laws
- The power handle will have longer life span if the above instructions are followed carefully

INSTRUMENT DESCRIPTION

The Acura® *electro* is a motorized, microprocessor controlled air displacement pipette. Energy is provided by a long life NiMH battery pack located in the pipette handle. This power handle can be replaced with a fully charged one within seconds, thus ensuring continuous work-flow without loss of set parameters. Acura® 925 models allows precise and reproducible pipetting within 0.25 µL to 1000 µL. 935 macromodels cover volumes from 0.25 to 10 mL. Multichannel Acura® 955 with 8 and 12 channel range between 0.5 µL and 350 µL. All volumetric modules (lower assemblies) are interchangeable on the same control unit (upper assembly). The Justip™ ejection system provides height adjustment of the single and multichannel tip ejector.

View (see p. I)

- A) Acura® *electro* pipettes
- B) Charging stand for electronic pipette and power handle
- C) Power supply with cord
- D) Power handle
- E) Charging LED
- F) Charging stand for power handle only (4.8V)

Keys and functions (fig. 1)

- 1) 2 positions start button
- 2) 3 positions speed selector
- 3) Programming key (mode) for choosing:
 - Pipetting modes
 - Side of display reading
 - Installed volumetric module
- 4) Set key for entering selection or emergency stop
- 5) Selection keys (+/-) for choosing:
 - Volumes
 - Calibration settings
 - Volume ranges
 - Left/right display reading
 - Mixing
- 6) LCD display (see fig. 3 for details)
- 7) Tip ejector button
- 8) Power handle
- 9) Control unit
- 10) Interchangeable volumetric module
- 11) Justip™ ejection system
- 12) Clips
- 13) Interchangeable multichannel module
- 14) Pipette nozzle

Start button (fig. 2A)

The start button has two contact positions:

- Press until first stop (1) to work at the slow pipetting speed
- Press until second stop (2) to work at the selected pipetting speed

LCD display (fig. 3)

- 1) Forward pipetting mode
- 2) Reverse pipetting mode
- 3) Dilution mode, volume sequence (vol. 1, 2 or 3)
- 4) Stepper mode
- 5) Low battery indicator
- 6) Pipetting steps: aspirating or dispensing
- 7) Digits for volume display or messages
- 8) Current volume units (µL or mL)
- 9) Sign requiring user input (selection or validation)

Power handle (fig. 4)

The rechargeable 4.8V power handle contains NiMH batteries of 300 mAh/4.8V. Charge a fully discharged power handle in less than 1.5 hours.

Optional charging stands

- Charging stand for 4.8V Acura® *electro* and power handles (fig. 6) with 3 charging positions
- Compact charging stand for 4.8V power handles only (fig. 5): allows to simultaneously charge up to 3 power handles

Power supply (fig. 7-8)

- Input: 100-240, 50/60 Hz
- Out put: 7.5 VDC
- Supplied with electrical cord and plug
- Various plug styles available depending on user location. See ordering information

USE OF THE INSTRUMENT

Supply contents

Accessories supplied with the Acura® *electro* may vary, depending on country. Check exact contents on packaging label. All elements also available separately. Refer to chapter “Ordering information” for more details.

Note: *Keep original packaging for adequate instrument protection during future transport or shipment.*

Inserting power handle (fig. 4A)

The 4.8V power handle fits at the back of the pipette control unit. Insert and click in as shown. The initialisation starts automatically, followed by a self-calibration test. The display shows “RE-CAL”.

Notes: *A 6V power handle cannot be mounted on a Acura electro 4.8V instrument.
Charge power handle completely before first use.*

Replacing power handle (fig. 4B)

To remove the handle from control unit, press lock trigger located at the bottom of the handle and lift gently.

Note: *Pipetting data are kept in memory even during battery change. They will show again automatically upon re-inserting power handle.*

Charging power handle (fig. 5 and 6)

The power handle can be charged in three different ways:

- 1) Attached to the *electro* pipette and placed on the 4.8V charging stand
- 2) Alone, placed on 4.8V charging stand of the *electro* pipette
- 3) Alone, placed on the 4.8V compact charging stand (to be ordered separately)

The red LED goes on when charging a power handle.

The green LED indicates that charging is completed and goes on stand by with minimal energy consumption.

Maximal battery capacity is obtained after a few full charging/discharging cycles.

At maximum capacity, power handles allow for over 3000 forward pipetting movements (full plunger stroke on a single channel micropipette) without recharging.











Notes: *It is not possible to charge a 4.8V power handle with a 6V charging stand (320.903 or 320.912). ONLY use charging stand 320.903.48 or 320.913.48.
If unused for 10 minutes, the Acura® *electro* switches automatically to an energy saving stand-by mode (display switches off) to guarantee a longer operation. Simply press start button to re-activate instrument.*

PROGRAMMING

Holding the electronic pipette (fig. 1)

The ergonomic shape of the Acura® *electro* allows long pipetting series without hand fatigue. Place finger rest on the phalanx of forefinger. The thumb reaches naturally start button (1) and ejector button (7), both easily activated. Multichannel casing (13) revolves to allow the selection of the best working position.

Left or right display reading










| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|---|--|---|--|
| Configuration for left or right handed operation | Press (> 0.5 sec.) |  |  |
| | Validate function |  |  |
| | Select |  |  |
| | Validate side |  |  |
| | Press (> 0.5 sec.) Back to current pipetting mode |  |  |

Note: After selection of left or right display reading, the last pipetting mode used will be displayed.

PROGRAMMING







Forward mode, programming and pipetting

In forward mode, the exact volume desired is aspirated. Liquid dispensing is automatically followed by a short excess plunger movement (blow out). Plunger returns back to its initial position one second later.

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--------------------------|-------------------|---|--|
| Forward pipetting | Select mode |  |  |
| | Validate mode |  |  <i>Last settings appear by default</i> |
| Choice of volume | Ex: select 400 μL |   |  |
| | Validate volume |  |  |

Pipetting in forward mode

Press start button gently until first stop to work at the slow pipetting speed. Press fully (second stop) to work at the selected pipetting speed (fig. 2A).










| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|-------------------|--------|---|--|
| Aspiration | Press |  |  |
| | | |  |
| Dispensing | Press |  |  |
| | | |  |

Notes: Plunger stays down when start button is kept pressed, and returns to initial position when released. Slightly touch reservoir wall when dispensing.

PROGRAMMING









Reverse mode, programming and pipetting

In reverse mode, the liquid is aspirated in excess of the selected volume. However, the set volume only will be delivered. The excess volume remains in the tip and can be kept or discarded. The reverse pipetting mode is recommended for viscous, volatile or foaming liquids.

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|---|-------------------|---|--|
| Choice of reverse pipetting mode | Select |  |  |
| | Validate mode |  |  <i>Last settings appear by default</i> |
| Choice of volume | Ex: select 150 μL |   |  |
| | Validate volume |  |  |

Pipetting in reverse mode

Press start button gently until first stop to work at the slow pipetting speed. Press fully (second stop) to work at the selected pipetting speed (fig. 2A).


| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|----------------------|--------------|---|--|
| Aspiration | Press |  |  |
| | | |  |
| Dispensing | Press |  |  |
| Excess volume | Double click |  |  |
| | | |  |

Notes: Skip "purge" by holding start button down after dosing is completed. The next sample is aspirated directly after release of the start button.
Slightly touch reservoir wall when dispensing.

PROGRAMMING

Stepper mode, programming and pipetting

In stepper mode, the volume aspirated is in small excess of the sum of every single aliquot. It is distributed step by step according to volume and number of aliquots programmed by the user.

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|---|-------------------|---|--|
| Choice of stepper pipetting mode | Select |  |  |
| | Validate mode |  |  |
| <i>Last settings appear by default</i> | | | |
| Choice of volume | Ex : select 50 μL |   |  |
| | Validate volume |  |  |
| <i>Display show maximal number of aliquots possible for the selected volume</i> | | | |
| Choice: number of aliquots | Ex: select 15 x |   |  |
| | Validate aliquots |  |  |






Pipetting in stepper mode

Press start button gently until first stop to work at the slow pipetting speed. Press fully (second stop) to work at selected pipetting speed (fig. 2A).

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|---------------------------------------|--------|---|--|
| Aspiration (Ex: 15 x 50 μL) | Press |  |  |
| | Ready | |  |
| Dispensing | Press |  |  |








PROGRAMMING

Pipetting in stepper mode (continued)

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|-------------------------------|-----------|---|---|
| Dispensing (continued) | | |  <p>Number of aliquots decrease after each dispensing</p> |
| Dispensing | Press |  |   |
| Excess volume | see below | |  |

Notes: Slightly touch reservoir wall when dispensing.

Stepper mode, excess volume control

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|---|--------------|---|---|
| Current display | | |  |
| Keep excess volume, ex. aspiration of same liquid | Press |  |   |
| or | | | |
| Blow out excess volume | Double click |  |  |
| Ready for aspiration of new liquid | | |  |

PROGRAMMING






















Stepper mode, interruption of pipetting sequence

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|-----------------------|--------------|------------|-------------|
| Current display | | | |
| End of pipetting | Press | STOP | |
| Residual volume | Double click | | |
| Ready for new filling | | | |

PROGRAMMING

Dilution mode, programming and pipetting

In dilution mode, up to 3 volumes of different liquids – each separated in the tip by an air bubble – will be aspirated. The total resulting volume is then dispensed in a single dose.







| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|-------------------|---|---|
| Choice of dilution pipetting mode | Select |  |  |
| | Validate mode |  |  <i>Last settings appear by default</i> |
| Choice of 1st volume | Ex: select 250 μL |   |  |
| | Validate volume 1 |  |  <i>Display shows volume capacity left to reach nominal volume</i> |
| Choice of 2nd volume | Ex: select 100 μL |   |  |
| | Validate volume 2 |  |  |
| Choice of 3rd volume | Ex : select 50 μL |   |  |
| No 3rd volume? | Select 0 μL |  |  |
| | Validate volume 3 |  |  |

Notes: When 3rd liquid aspirated (V3), air bubble prohibits reaching cumulated volume equal to maximum pipette capacity. Thus avoiding risk of nozzle contamination.
Air bubbles in macrotips (Acura® 935 models) only serve the purpose of separating liquid at the lower tip area.

PROGRAMMING

Pipetting in dilution mode

Press start button gently until first stop to work at the slow pipetting speed. Press fully (second stop) to work at the selected pipetting speed (fig. 2A).

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|--------------------------------------|---|--|
| Aspiration 1st volume | Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL V 1</small> ASP <small>µl</small> </div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center; margin-top: 5px;"> <small>DIL</small> AIR </div> |
| Aspiration air bubble | Lift tip out of the liquid, Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL</small> ASP </div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center; margin-top: 5px;"> <small>DIL V 2</small> <small>ASP</small> 100 <small>µl</small> </div> |
| Aspiration 2nd volume | Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL V 2</small> <small>ASP</small> ASP <small>µl</small> </div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center; margin-top: 5px;"> <small>DIL</small> AIR </div> |
| Aspiration air bubble | Lift tip out of the liquid, Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL</small> ASP </div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center; margin-top: 5px;"> <small>DIL V 3</small> <small>ASP</small> 50 <small>µl</small> </div> |
| Aspiration 3rd volume <i>(if programmed)</i> | Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL V 3</small> <small>ASP</small> ASP <small>µl</small> </div> <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center; margin-top: 5px;"> <small>DISP</small> <small>DIL</small> 400 <small>µl</small> </div> |
| Dispensing V1+V2+V3 | Press |  | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL</small> <small>DISP</small> DISP <small>µl</small> </div> |
| Ready for new filling | | | <div style="border: 1px solid black; border-radius: 10px; padding: 5px; text-align: center;"> <small>DIL V 1</small> <small>ASP</small> 250 <small>µl</small> </div> |

Note: Slightly touch reservoir wall when dispensing.

PROGRAMMING






Dilution mode, interruption of pipetting sequence

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--------------------------|--------------|---|--|
| Current display | | |  |
| End of pipetting | Press |  |  |
| Blow out residual volume | Double click |  |  |
| Ready for new filling | | |  |

PROGRAMMING

Mixing

Available in all pipetting modes, the mixing performs consecutive back and forth aspiration/dispensing of the set volume. It is only possible after pipetting or purge steps are over.

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|------------------|---|--|
| Mixing with 3 cycles of aspiration/dispensing | Press 1 x |  |  |
| Continuous mixing | Keep key pressed |  |  |
| Back to the pipetting | Press |  |  |




OPERATION

Selection of pipetting speed (fig. 2B)

Three pipetting speeds are available. Simply move the selector from left to right (slow/fast) to adapt working speed to the type of liquid or to a specific application. Independently from the speed selection, it is always possible to aspirate or distribute at the slow speed by gently pressing start button half way (feel first stop).

Charge level of power handle (fig. 3)

Pay attention to charge indicator on LCD to avoid unexpected power cut-off. For appropriate handling of the power handle, refer to chapters on “Replacing power handle” (fig. 4B) and “Charging power handle” (fig. 5 and 6) for more details.

| LCD DISPLAY | SIGNIFICATION/CAUSE | SOLUTION |
|---|----------------------|--|
|  | Charged battery | Work without trouble |
|  | Low battery | It is recommended to charge battery after pipetting work is finished or to replace power handle |
|  | Battery almost empty | Energy level insufficient Instrument will switch off Replace or charge power handle immediately with a 4.8V power handle |

Note: Always keep one or more spare 4.8V power handle(s) available on the charging stand.

Changing volumetric module

Lower assembly of the Acura® *electro* is interchangeable on a single control unit (upper assembly). Pipetting sequence must be completed before disassembling.

Note: 350 μL volumetric modules are only compatible with Acura® *electro* pipettes bearing serial numbers xxxx3001 or higher.

Disassembling volumetric module up to 1000 μL (fig. 9 inside front cover)

- ① Press ejector button to the bottom
- ② Hold screw of tip ejector and slightly turn to the left. Pull out tip ejector (see arrows on handle).
- ③ Unscrew barrel
- ④ Unclip plunger

Before storing volumetric module, reassemble plunger, barrel, and ejector and place protection cap (accessory, Cat. No. 825.691) on the top.

Disassembling volumetric module 5 and 10 mL (fig. 10 inside front cover)

- ① Slightly turn ejector nut to unclip from ejector cap
- ② Unscrew barrel, then gently pull barrel to unclip plunger rod
- ③ Press ejector button to the bottom
- ④ Turn ejector cap to the left and pull out

Before storing volumetric module, reassemble ejector nut and cap then place protection cap (accessory, Cat. No. 825.691) on the top.

Assembling

Remove protection cap and pull out plunger

- ⑤ Hold plunger between thumb and index finger and clip plunger

Note: *Plunger must be clicked in the control unit before assembling the module.*

- ③ Screw on barrel

- ① Press ejector button to the bottom

- ② Introduce ejector screw into the control unit (see arrows on handle). Slightly turn right to lock. Release ejector button

Enter module data in control unit prior to use, as described in chapter “Programming volumetric module”.

Assembling

Remove protection cap

- ⑤ Pull out plunger rod and introduce tip or rod in side hole to prevent from retracting
- ⑥ Introduce plunger rod in control unit and click in

Note: *Plunger must be clicked in the control unit before assembling the module. If difficulty to pull out plunger, use the small pin supplied in the box when ordering additional volumetric module.*

Hold pipette nozzle and screw barrel firmly

- ③ Press ejector button, introduce tooth of ejector cap in recess (see arrows on handle), and turn right to lock

Enter module data in control unit prior to use as described in chapter “Programming volumetric module”.

Changing volumetric module *(continued)*

Disassembling multichannel volumetric module (fig. 12 inside front cover)

- ① Press ejector button
- ② Turn ejector nut to the left till its lowest position, unclip ejector ring and release ejector button
- ③ Hold volumetric module firmly with one hand and, while pulling down casing, rotate slowly to unscrew
- ④ Gently unclip plunger rod.

Before storing volumetric module, place protection cap (accessory, Cat. No. 825.691) on the top.

Assembling

Remove protection cap

- ⑤ Pull out plunger rod and introduce small pin (i.e. paper clip) in side hole to prevent from retracting
- ⑥ Introduce plunger rod in control unit and click in

Note: Plunger must be clicked in the control unit before assembling the module.

Hold casing with fingers pulling against barrel extremities

Slowly screw volumetric module while positioning tooth of ejector ring in recess (see arrows on handle)

- ③ Press ejector button, turn ejector ring to the right until tooth clicks in and release ejector button

Control proper functioning of ejector button. Set ejector height to desired position. Enter module data in control unit prior to use as described in chapter "Programming volumetric module".

Notes: Power handle **must** be connected on control unit prior to re-assembling volumetric module. Watch for not activating programming keys when changing volumetric module.

IMPORTANT, before the very first use of a new volumetric module, perform pipette calibration according to chapter "Calibration".

OPERATION

Programming for volumetric module

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|--------|------------|-------------|
| Choice of volumetric module | | | 1000? |
| <i>Last settings appear by default</i> | | | |
| Ex: module ranging 20-200 μL | | + - | 200? |
| Validate the choice | | SET | RE-CAL |
| <i>Instrument performs a self-calibration test automatically</i> | | | |

| | | | |
|--|--|--|----------|
| Select pipetting mode according to chapter "Programming" | | | FORWARD? |
|--|--|--|----------|




Correction of volumetric module selection

It is also possible to access the menu "Programming volumetric module" at anytime.

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|---|------------|------------------------------|
| Current display | | | PIP ASP 100 μL |
| Choice of volumetric module menu | Press (> 0.5 sec.) | MODE | SIDE ? |
| | Press | MODE | MODULE? |
| | Validate | SET | 200? μL |
| Correction of volumetric module selection ▼ | Ex: 8-channel module ranging 5-50 μL | + - | 8x50? μL |

OPERATION

Correction of volumetric module selection *(continued)*

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|--------|---|--|
| Correction of volumetric module <i>(continued)</i> | |  |  <i>Instrument performs self-calibration test automatically</i> |
| Select pipetting mode according to chapter “Programming” | | |  |

Note: WARNING, the programmed volumetric module must correspond to the one fitted on the pipette assembly.

Use of Pasteur pipette (935, 5 mL model only)

Glass Pasteur pipettes are of advantage in handling PP affecting solvents. Optional adapter nozzle (cat. No. 1.835.633) fits 0.5-5 mL model to accommodate standard 2 mL Pasteur pipettes (ext. Ø 6.5-7.2 mm) in addition to Polypropylene tips. Refer to chapter “Replacement of tightness parts” for disassembling and changing adapter nozzle. Slightly grease O-rings in adapter for tight Pasteur pipette fitting.



Note: WARNING, do not set volume higher than 2 mL.

Adjustment of tip ejector (fig. 14 and 15)

The Justip™ system allows instant height adjustment (+/- 2 mm) of ejector for best fitting of the tip used (optimal distance).

Set the correct position (~0.5 mm space between tip and ejector) by rotating ejector screw to the left or to the right (LO-HI). Click stops prevent any unwanted change while pipetting.

Note: Ejector head on multichannel models designed for soft, sequential tip removal.



MAINTENANCE AND STERILISATION

The Acura® *electro* requires little basic maintenance. To obtain long term trouble free use, special attention during pipetting and regular cleaning are recommended.

Cleaning

- External parts of control unit, power handle and charging stands are cleaned with dampened cloth
- Volumetric module, once disassembled according to instructions in chapter “Operation”, can be cleaned with alcohol, or soaked in decontamination or disinfecting solutions. Ultra-sonic bath helps to remove sticking residues
- On volumetric modules up to 1000 µL, plunger tightness is provided by PTFE sleeve on O-ring. Grease O-ring slightly after sleeve was disassembled. Also grease plunger, O-ring and barrel wall of 5 and 10 mL models before reassembling
- Any defective part must be replaced. Order original spare parts from authorised dealers

Note: *WARNING: No liquid must penetrate into control unit (upper assembly).*

Replacement of tightness parts

O-ring and PTFE sleeve, single channel micro-volumetric modules (up to 1000 µL)

- Remove volumetric module from control unit according to instructions in chapter “Operation” (fig. 9)
- Press both clips of the ring (Cat. No. 825.701) with a pipette tip or a pointed tool. (See exploded drawings p. IV)
- Remove ring from barrel
- Remove spring to access O-ring/PTFE sleeve assembly
- Change parts. Clean plunger to eliminate any grease residue. Slightly grease O-ring and reassemble according to instructions in chapter “Operation”

Note: *Tightness parts are not available separately on 10 µL model. In case of tightness problem, barrel must be changed. Plunger must be “clicked” with the control unit **before** assembling the volumetric volume.*

O-ring, single channel macro-volumetric modules (5 and 10 mL)

- Remove volumetric module from control unit according to instructions in chapter “Operation” (fig. 10)
- Press both barrel clips with fingers to separate from bonnet (fig. 11)
- Pull out plunger assembly. Unscrew plunger rod and remove washers and spring
- Change parts if needed. Grease O-ring, washer and barrel
- Reassemble plunger assembly, barrel and bonnet according to instructions in chapter “Operation”

Note: *Plunger must be “clicked” with the control unit **before** assembling the volumetric volume.*

Replacement of tightness parts *(continued)*

Barrel change, multichannel volumetric modules (up to 350 µL)

Note: O-ring cannot be removed from barrel. Change barrel if tightness is deficient.

- Remove volumetric module from control unit according to instructions in chapter «Operation» (fig. 12)
- Press both clips of cover with pointed tool and remove casing (fig. 13)
- Press bottom clips (a) of barrel holder and separate (b)
- Pull out barrels (c)
- Apply thin, even grease layer on all surface length of plunger before reassembling
- Introduce barrel on plunger. Separate both plate of barrel holder, reposition barrel
- Clip barrel plate, all barrels must be sitting properly and aligned
- Place assembly in casing and clip on cover

Notes: Plunger must be "clicked" with the control unit **before** assembling the volumetric volume. Markings on casing (volume) and on cover (Justip) should appear on opposite sides.

Sterilisation

Only the volumetric module (lower part) is autoclavable at 121°C (20 minutes, 1 atm). Disassemble from control unit according to instructions in chapter "Operation". Before autoclaving, remove nozzle filter on 935 models. Autoclave volumetric module as one assembly. Parts must be cooled down and completely dry before mounting on control unit. Check tightness and accuracy after first cycle, then regularly but at least after 50 autoclaving cycles. Correct autoclaving and resulting sterility are the responsibility of the user.

Note: Setting parameters on the Acura® electro must correspond to those of the assembled volumetric module.

INSTRUMENT CALIBRATION

Each Acura® *electro* pipette is factory tested for conformity according to ISO 8655 standards. Calibration parameters are permanently memorised in the instrument's microprocessor. If performance results are no longer within recommended values, for instance after QC check, or replacing parts, or if changes occur in physical parameters (liquid density, temperature, atmospheric pressure), re-calibration is easily performed through the calibration menu.







Notes: *WARNING, a calibration is necessary before the very first use of a volumetric module other than the one supplied with the instrument (even if same range). It is recommended to control instrument performance in accordance with internal laboratory procedures (SOP/GLP, etc.) or at least once a year.*

Calibration increments

In the calibration menu, "QC-CAL" is the original calibration value of a factory-calibrated instrument. Other volumetric modules purchased as accessories will display "0" as target calibration value. Calibration is performed by changing the calibration unit figure according to following values:

| | | | | | | | |
|------------------------------------|-------------|-------------|-----------|----------|-----------|-----------|----------|
| Volumetric module (lower assembly) | 10 µL | 50 µL | 200 µL | 350 µL | 1000 µL | 5 mL | 10 mL |
| Each calibration increment = | ± 0.0025 µL | ± 0.0125 µL | ± 0.05 µL | ± 0.1 µL | ± 0.25 µL | ± 1.25 µL | ± 2.5 µL |

New calibration

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|----------------------------------|---------------------------|---|--|
| Finish pipetting sequence | | |  |
| Access calibration menu | Press first |  | |
| | Then press simultaneously |  |  |
| | Validate |  |  |

Display shows "QC-CAL" for a factory-calibrated instrument where the original parameters have not been modified. If modified, display shows the number of calibration increments previously selected


















INSTRUMENT CALIBRATION

| New calibration (continued) | | | |
|------------------------------|--|------------|------------------------|
| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
| Change calibration parameter | Ex: volume reduction of 0.75 μL (= 3 x 0.25 μL / increment) on a 1000 μL model Press | + - | - 3? |
| | Validate | SET | - 3 |
| | Press | MODE | PIP ASP 4 0 0 μ |

Note: When purchasing a new volumetric module as an accessory, it is important to introduce calibration parameters before very first use. To do so, perform calibration according to chapter "Calibration". The value entered is automatically and permanently memorised in the microprocessor when leaving calibration menu.




TROUBLE SHOOTING

| Error messages | | |
|---|---|---|
| LCD DISPLAY | CAUSE | SOLUTION |
|  | Plunger sticking or dragging | Disassemble volumetric module according to chapter "Operation" Clean plunger according to chapter "Maintenance" Reset instrument (see below) |
| | Plunger rod not clipped before assembling | Reconnect plunger rod according step 5 and 6 of page 18 and 19. Reset instrument (see below) Contact authorised dealer for control if error persists |
|  | Microprocessor detected deviation between set volume and effective plunger travel | Reset instrument (see below) |
| | Plunger rod not clipped before assembling | Reconnect plunger rod according step 5 and 6 of page 18 and 19. Reset instrument (see below) |
|  | Volumetric module disassembled when pipetting sequence is not completed | Reset instrument (see below) Confirm selection of volumetric module |
|  | Only in dilution mode Volumes programmed larger than maximal aspiration capacity | Reset instrument (see below) |

| Instrument reset | | | |
|-------------------|------------------------------|---|--|
| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
| Error message | | |  |
| | Press | STOP  |  |
| Reset | Double click on start button |  |  |
| Error message | | |  |
| | Press | STOP  |  |
| Select new volume | Press |   |  |

TROUBLE SHOOTING

Instrument reset *(continued)*

| OPERATION | ACTION | KEY/BUTTON | LCD DISPLAY |
|--|------------|---|--|
| Or change calibration settings | Press then |   |  |
| Refer to chapter "Calibration" for new setting | | | |

Other failures

| Observation | Possible cause | Action |
|---|---|---|
| Power handle does not fit in pipette control unit | Control if compatible 4.8V device | Exchange with 4.8V power handle |
| Power handle or instrument is not charging/no red light | Control if compatible 4.8V device | Exchange with 4.8V device |
| Plug do not fit into stand | Control for correct power supply model 4.8V | Use 4.8V power supply |
| No display | Instrument in stand-by | Press start button to activate instrument |
| | Battery is discharged | Charge power handle or replace if damaged |
| LCD display on but no reaction when pressing start button | Volumetric module not correctly locked | Check volumetric module |
| Poor instrument performance | Lack of tightness | Check proper tip fitting. Use tips compatible with the instrument Check tip cone, change if damaged Check O-ring and PTFE sleeve, change if damaged |
| | Instrument not calibrated | Perform new calibration |
| | Instrument pipetting viscous, or volatile solutions, liquid temperature not comprised between 20-25°C | New calibration with specific solution or temperature |
| Reduced battery life span | Power handle damaged | Replace 4.8V power handle |
| | Too much friction in lower assembly | Clean lower assembly |
| Red light on charger stand not lid | Micropipette or power handle not sitting correctly on stand | Reposition instrument or power handle |
| | Using a 6V power handle instead of 4.8V | Replace with a 4.8V power handle |
| Wrong volume | Erroneous programming of volumetric module | Set parameters correctly |
| Long pipetting time Plunger moves stepwise | Plunger sticking or dragging | Disassemble volumetric module and clean |
| | Motor drive impaired | Contact authorized dealer for control |

PERFORMANCE

Accepted tolerance value as obtained with bi-dist. water at constant temperature ($\pm 0.5^{\circ}\text{C}$) comprised between 20°C and 25°C , according to ISO 8655. Use original or compatible tips.

Acura® electro 925

| Volume μL | Division μL | Tip type μL | Performance at μL | Inaccuracy E% | Imprecision CV% |
|-------------------------|---------------------------|---------------------------|---------------------------------|------------------|--------------------|
| 0.5 - 10 | 0.05 | 10 | 1 | $< \pm 2.5 \%$ | $< 1.8 \%$ |
| | | | 5 | $< \pm 1.8 \%$ | $< 1.2 \%$ |
| | | | 10 | $< \pm 1.0 \%$ | $< 0.5 \%$ |
| 2.5 - 50 | 0.25 | 200 | 5 | $< \pm 2.5 \%$ | $< 1.5 \%$ |
| | | | 25 | $< \pm 1.6 \%$ | $< 0.9 \%$ |
| | | | 50 | $< \pm 0.7 \%$ | $< 0.3 \%$ |
| 10 - 200 | 1.00 | 200 | 20 | $< \pm 1.5 \%$ | $< 0.7 \%$ |
| | | | 100 | $< \pm 1.1 \%$ | $< 0.5 \%$ |
| | | | 200 | $< \pm 0.6 \%$ | $< 0.2 \%$ |
| 50 - 1000 | 5.00 | 1000 | 100 | $< \pm 1.5 \%$ | $< 0.5 \%$ |
| | | | 500 | $< \pm 1.0 \%$ | $< 0.4 \%$ |
| | | | 1000 | $< \pm 0.5 \%$ | $< 0.2 \%$ |

Acura® electro 935

| Volume mL | Division mL | Tip type mL | Performance at mL | Inaccuracy* E% | Imprecision* CV% |
|--------------|----------------|----------------|----------------------|-------------------|---------------------|
| 0.25 - 5 | 0.05 | 5 | 0.5 | $< \pm 1.0 \%$ | $< 0.6 \%$ |
| | | | 2.5 | $< \pm 0.9 \%$ | $< 0.5 \%$ |
| | | | 5 | $< \pm 0.7 \%$ | $< 0.3 \%$ |
| 0.5 - 10 | 0.05 | 10 | 1 | $< \pm 0.8 \%$ | $< 0.4 \%$ |
| | | | 5 | $< \pm 0.7 \%$ | $< 0.3 \%$ |
| | | | 10 | $< \pm 0.5 \%$ | $< 0.2 \%$ |

* Macropipette with nozzle filter.

Acura® electro 955 (8 and 12-channel)

| Volume μL | Division μL | Tip type μL | Performance at μL | Inaccuracy E% | Imprecision CV% |
|-------------------------|---------------------------|---------------------------|---------------------------------|------------------|--------------------|
| 0.5 - 10 | 0.05 | 10 | 1 | $< \pm 3.5 \%$ | $< 3.0 \%$ |
| | | | 5 | $< \pm 2.5 \%$ | $< 2.0 \%$ |
| | | | 10 | $< \pm 1.5 \%$ | $< 1.0 \%$ |
| 2.5 - 50 | 0.25 | 200 | 5 | $< \pm 1.0 \%$ | $< 1.0 \%$ |
| | | | 25 | $< \pm 0.9 \%$ | $< 0.7 \%$ |
| | | | 50 | $< \pm 0.8 \%$ | $< 0.4 \%$ |
| 20 - 350 | 5.00 | 350 | 40 | $< \pm 1.0 \%$ | $< 0.6 \%$ |
| | | | 200 | $< \pm 0.9 \%$ | $< 0.4 \%$ |
| | | | 350 | $< \pm 0.8 \%$ | $< 0.3 \%$ |

Notes: Use of other tips than those recommended, as well as pipetting viscous or volatile liquids may lead to performance deviation compared to those shown in the above figure. Product specifications subject to change without prior notice.

Performance value obtained in forward mode. Small deviations may exist when using other pipetting modes. To obtain best possible performance with one specific pipetting mode, it is recommended to perform a new calibration.

WARRANTY

Your Acura® *electro* and power handle are guaranteed against any material or manufacturing defects for the period of time specified in its QC certificate. Damages due to non-respect of manufacturer's instructions, safety precautions or autoclaving conditions, as well as material colour alteration are excluded from the warranty. Repair and replacement of parts do not extend warranty time. Claims for warranty are void if instrument has been tempered. Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased after obtaining return authorisation.

Note: Decontaminate volumetric module of the instrument prior to returning it.



ORDERING INFORMATION

Micropipettes

Acura® electro 925

| Volume µL | Division µL | Tip type µL | Pipette alone Cat. No. | Basic package* Cat. No.** |
|--------------|----------------|----------------|---------------------------|------------------------------|
| 0.5 – 10 | 0.05 | 10 | 925.0010.48 | 925.0010.48x |
| 2.5 – 50 | 0.25 | 200 | 925.0050.48 | 925.0050.48x |
| 10 – 200 | 1.00 | 200 | 925.0200.48 | 925.0200.48x |
| 50 – 1000 | 5.00 | 1000 | 925.1000.48 | 925.1000.48x |

Acura® electro 935

| Volume mL | Division mL | Tip type mL | Pipette alone Cat. No. | Basic package* Cat. No.** |
|--------------|----------------|----------------|---------------------------|------------------------------|
| 0.25 – 5 | 0.05 | 5 | 935.05.48 | 935.05.48x |
| 0.5 – 10 | 0.05 | 10 | 935.10.48 | 935.10.48x |

Acura® electro 955

| Volume µL | Division µL | No of channel | Tip type µL | Pipette alone Cat. No. | Basic package* Cat. No.** |
|--------------|----------------|------------------|----------------|---------------------------|------------------------------|
| 0.5 – 10 | 0.05 | 8 | 10 | 955.08.010.48 | 955.08.010.48x |
| 2.5 – 50 | 0.25 | 8 | 200 | 955.08.050.48 | 955.08.050.48x |
| 20 – 350 | 5.00 | 8 | 350 | 955.08.350.48 | 955.08.350.48x |
| 0.5 – 10 | 0.05 | 12 | 10 | 955.12.010.48 | 955.12.010.48x |
| 2.5 – 50 | 0.25 | 12 | 200 | 955.12.050.48 | 955.12.050.48x |
| 20 – 350 | 5.00 | 12 | 350 | 955.12.350.48 | 955.12.350.48x |

* Basic package includes electronic micropipette with individual QC certificate, pipette charging stand, power supply with cord, an additional power handle, Qualitips pipette tips samples and operating instructions.

** Replace x by letter for correct plug type depending on user location: E=Europe, G=UK, U=USA-Japan, A=Australia/NZ.

Accessories and stands

| Description | Packaging | Cat. No. |
|---|-----------|------------|
| Interchangeable, autoclavable single channel volumetric module | | |
| 0.5 – 10 µL | 1/pk | 800.0010 |
| 2.5 – 50 µL | 1/pk | 800.0050 |
| 10 – 200 µL | 1/pk | 800.0200 |
| 50 – 1000 µL | 1/pk | 800.1000 |
| 0.25 – 5 mL | 1/pk | 800.5000 |
| 0.5 – 10 mL | 1/pk | 800.10000 |
| Interchangeable, autoclavable 8-channel volumetric module | | |
| 0.5 – 10 µL | 1/pk | 800.08.010 |
| 2.5 – 50 µL | 1/pk | 800.08.050 |
| 20 – 350 µL | 1/pk | 800.08.350 |
| Interchangeable, autoclavable 12-channel volumetric module | | |
| 0.5 – 10 µL | 1/pk | 800.12.010 |
| 2.5 – 50 µL | 1/pk | 800.12.050 |
| 20 – 350 µL | 1/pk | 800.12.350 |

ORDERING INFORMATION

Accessories and stands *(continued)*

| Description | Packaging | Cat. No. |
|--|-----------|-------------|
| Power handle (battery pack) | 1/pk | 900.910.48 |
| Power handle (duo battery pack) | 2/pk | 900.912.48 |
| Charging stand for <i>electro</i> pipettes and power handles* | 1/pk | 320.903.48 |
| Compact charging stand for power handles only* | 1/pk | 320.912.48 |
| <i>* Power supply with plug to be ordered separately</i> | | |
| Power supply 100-240 V with cord and plug for Europe | 1/pk | 900.901.48E |
| Power supply 100-240 V with cord and plug for UK | 1/pk | 900.901.48G |
| Power supply 100-240 V with cord and plug for USA-Japan | 1/pk | 900.901.48U |
| Power supply 100-240 V with cord and plug for Australia/NZ | 1/pk | 900.901.48A |
| Universal work station 337 for 7 single channel pipettes or spare volumetric modules | | |
| Light grey | 1/pk | 320.337G |
| Ice blue | 1/pk | 320.337B |
| Vanilla yellow | 1/pk | 320.337Y |
| Mint green | 1/pk | 320.337M |
| Pastel rose | 1/pk | 320.337R |

Qualitips® pipette tips

| Description | Packaging | Cat. No. |
|---|---------------|------------|
| Selection of Qualitips® without filter | | |
| Ultra-micro tips, 10 µL | 1x1000/bag | 309.0010B |
| | 10x 96/rack | 309.0010R |
| Natural tips, 200 µL | 2x 500/bag | 307.0200B |
| | 10x 96/rack | 307.0200R |
| Yellow tips, 200 µL | 2x 500/bag | 327.0200B |
| | 10x 96/rack | 327.0200R |
| Natural tips, 200 µL | 1x1000/bag | 308.0200B |
| | 10x 96/rack | 308.0200R |
| Natural tips, 350 µL | 1x1000/bag | 308.0350B |
| | 10x 96/rack | 308.0350R |
| Blue tips, 1000 µL | 4x 250/bag | 318.1000B |
| | 8x 96/rack | 318.1000R |
| Natural tips, 5 mL | 2x 250/bag | 312.05B |
| | 2x 50/rack | 312.05R |
| Natural tips, 10 mL | 2x 250/bag | 312.10 |
| | 4x empty rack | 312.10ER |
| Selection of Qualitips® with filter | | |
| Ultra-micro tips, 10 µL, sterile | 10x 96/rack | 309.0010FR |
| Natural tips, 100 µL, sterile | 10x 96/rack | 308.0100FR |
| Natural tips, 300 µL, sterile | 10x 96/rack | 308.0300FR |
| Natural tips, 1000 µL, sterile | 10x 100/rack | 309.1000FR |
| Nozzle protection filter, cellulose | | |
| For 5 mL volumetric module | 1x 250/bag | 322.05 |
| For 10 mL volumetric module | 1x 100/bag | 322.10 |




















Note: Refer to sales literature for full programme of Qualitips® pipette tips.

0.5 - 10 μ L

2.5 - 50 μ L
















10 - 200 μ L

50 - 1000 μ L



0.25 - 5 mL
















0.5 - 10 mL



0.5 - 10 μ L

2.5 - 50 μ L

20 - 350 μ L



Operating instructions also available in other languages



U. S. S.
F C C
part 15



ENERGY STAR®
qualified power
supply for a better
environment



SAS Accredited
Service Laboratory
ISO/IEC 17025



SOCOREX ISBA S.A. – CH. CHAMP-COLOMB 7 – 1024 ECUBLENS/LAUSANNE – SWITZERLAND
socorex@socorex.com – www.socorex.com – TEL. +41 21 651 6000 – FAX +41 21 651 6001