



Qualitix®

The new Qualitix® pipette tip programme is all about quality, flexibility and compatibility. An extended line made of 100% pure virgin, medical grade polypropylene, fitting not only Socorex® but all major micropipette brands. Packaging selection (bags, racks and *tipfill* rack filling system), together with a large selection of filter tips, offers appreciated versatility.

Free samples available on request.


The Qualitix® pipette tips from 10 µL to 10 mL

Product advantages

- Wide selection, superior quality, autoclavable
- Bags, racks and *tipfill* rack filling system
- Lot by lot consistency
- Certified free of nuclease, endotoxin and ATP
- Reduced packaging size, less bench space requirement
- Broad compatibility with other micropipette brands

Extended certifications

Every lot certified nuclease (DNase, RNase), endotoxin (pyrogen) and ATP-free. In addition, sterile tips guaranteed free of human DNA. Certificate supplied in each pack. CE conform.



Quick rack filling

tipfill system is intended for easy rack loading. Tip trays are instantly transferred into a rack without risk of contamination.

The smart re-closable *tipfill* system box reduces storage space, waste volume and shipping costs.

Low retention – maximal recovery




Low binding tips are made with a silicone-free resin reducing surface adhesion. They provide maximal liquid recovery. Best suited for proteomic research and handling of any viscous liquids. Available in 10, 200 and 1000 µL capacities, in sterile or non-sterile racks. Protease-free certified.



Aerosol protection – no contamination



Filter tips include a hydrophobic HDPE filter with 10 µ pore size to block any contaminant. Very high purity level certified.



Empty tip racks can be recycled and are compatible with *tipfill* system.





Ordering information – Qualitix® pipette tip selection

10 µL Regular Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	31.4 mm	Natural	10 µL ultra-microtips, fine pointed		1000 / bag	10 bags	3090.0010B
		Natural	10 µL ultra-microtips, fine pointed		12x 96 / rack	10 packs	3090.0010R
		Natural	10 µL <i>tipfill</i> system incl. ultra-microtips, fine pointed		13x 96 / pack	10 packs	3090.0010T
		Natural	10 µL ultra-microtips, fine pointed	✓	12x 96 / rack	10 packs	3090.0010RS
			Empty rack for 10/200 µL tips to be used with <i>tipfill</i> system		10 / pack	10 packs	3070.0210ER
	38.1 mm	Natural	10 µL extended microtips, fine pointed		1000 / bag	10 bags	3020.0010B
		Natural	10 µL extended microtips, fine pointed		10x 96 / rack	10 packs	3020.0010R
	71.0 mm	Natural	10 µL capillary tips for gel loading		4x 204 / rack	N/A	302.0010GR

10 µL Filter Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	31.4 mm	Natural	10 µL ultra-microtips, fine pointed	✓	12x 96 / rack	10 packs	3090.0010FRS
	38.1 mm	Natural	10 µL extended microtips, fine pointed	✓	10x 96 / rack	10 packs	3020.0010FRS



10 µL Low Binding Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	31.4 mm	Natural	10 µL ultra-microtips, low binding, fine pointed		12x 96 / rack	10 packs	3090.0010PR
		Natural	10 µL ultra-microtips, low binding, fine pointed	✓	12x 96 / rack	10 packs	3090.0010PRS
	38.1 mm	Natural	10 µL extended microtips, low binding, fine pointed		10x 96 / rack	10 packs	3020.0010PR
		Natural	10 µL extended microtips, low binding, fine pointed	✓	10x 96 / rack	10 packs	3020.0010PRS


20 - 100 µL Filter Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	49.6 mm	Natural	20 µL filter tips, bevel point	✓	10x 96 / rack	10 packs	3080.0020FRS
	49.6 mm	Natural	100 µL filter tips, bevel point	✓	10x 96 / rack	10 packs	3080.0100FRS


200 µL Regular Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	49.6 mm	Natural	200 µL universal tips, bevel point		1000 / bag	10 bags	3070.0200B
		Natural	200 µL universal tips, bevel point		10x 96 / rack	10 packs	3070.0200R
		Natural	200 µL <i>tipfill</i> system incl. universal tips, bevel point		10x 96 / pack	10 packs	3070.0200T
		Natural	200 µL universal tips, bevel point	✓	10x 96 / pack	10 packs	3070.0200RS
			Empty rack for 10/200 µL tips to be used with <i>tipfill</i> system		10 / pack	10 packs	3070.0210ER
	49.6 mm	Yellow	200 µL regular tips, bevel point		1000 / bag	10 bags	3290.0200B
		Yellow	200 µL regular tips, bevel point		10x 96 / rack	10 packs	3290.0200R
		Yellow	200 µL regular tips, bevel point	✓	10x 96 / rack	10 packs	3290.0200RS
	49.6 mm	Yellow	200 µL short neck tips for Calibra pipettes		2x 500 / bag	N/A	328.0200B
		Yellow	200 µL short neck tips for Calibra pipettes		10x 96 / rack	N/A	328.0200R
	82.5 mm	Natural	200 µL capillary tips for gel loading		6x 96 / rack	8 packs	3080.0200GR


Note: pack is the regular selling unit. When ordering the number of packs corresponding to a full case, Socorex supplies one case with an attractive rebate.





Ordering information – Qualitix® pipette tip selection


200 µL Filter Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	51.2 mm	Natural	200 µL universal tips, bevel point	✓	10x 96 / rack	10 packs	3070.0200FRS
	83.8 mm	Natural	200 µL extended tips, fine pointed	✓	6x 96 / rack	8 packs	3080.0200LFRS


200 µL Low Binding Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	49.6 mm	Natural	200 µL universal tips, low binding		10x 96 / rack	10 packs	3070.0200PR
		Natural	200 µL universal tips, low binding	✓	10x 96 / rack	10 packs	3070.0200PRS




350 µL Regular Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	58.8 mm	Natural	350 µL regular tips		1000 / bag	N/A	308.0350B
		Natural	350 µL regular tips		10x 96 / rack	N/A	308.0350R

300 µL Filter Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	58.8 mm	Natural	300 µL regular tips	✓	10x 96 / rack	N/A	308.0300FR

1000 µL Regular Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	88.8 mm	Natural	1000 µL universal tips, fine pointed		1000 / bag	10 bags	3070.1000B
		Natural	1000 µL universal tips, fine pointed		6x 96 / rack	8 packs	3070.1000R
		Natural	1000 µL tipfill system incl. universal tips, fine pointed		5x 96 / pack	10 packs	3070.1000T
		Natural	1000 µL universal tips, fine pointed	✓	6x 96 / rack	8 packs	3070.1000RS
			Empty rack for 1000 µL tips to be used with tipfill system		6 / pack	8 packs	3070.1000ER
	75.9 mm	Blue	1000 µL regular tips, bevel point		1000 / bag	10 bags	3190.1000B
		Blue	1000 µL regular tips, bevel point		6x 96 / rack	8 packs	3190.1000R
		Blue	1000 µL regular tips, bevel point	✓	6x 96 / rack	8 packs	3190.1000RS
	75.9 mm	Natural	1000 µL regular tips, bevel point		1000 / bag	10 bags	3090.1000B
		Natural	1000 µL regular tips, bevel point		6x 96 / rack	8 packs	3090.1000R
	72.0 mm	Blue	1000 µL regular tips		4x 250 / bag	N/A	319.1000B
		Blue	1000 µL regular tips		8x 60 / rack	N/A	319.1000R

1000 µL Filter Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	88.9 mm	Natural	1000 µL universal tips, fine pointed	✓	6x 96 / rack	8 packs	3070.1000FRS

1000 µL Low Binding Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	88.9 mm	Natural	1000 µL universal tips, low binding, fine pointed		6x 96 / rack	8 packs	3070.1000PR
		Natural	1000 µL universal tips, low binding, fine pointed	✓	6x 96 / rack	8 packs	3070.1000PRS

2, 5, 10 mL Regular Tips	Size	Description			Packaging		Cat. No.
		Colour	Style	Sterile	Pack	Case	
	115.5 mm	Natural	2 mL macrotips		2x 250 / bag	N/A	312.02
	123.2 mm	Natural	5 mL macrotips, graduated		2x 250 / bag	N/A	312.05B
		Natural	5 mL macrotips, graduated		2x 50 / rack	N/A	312.05R
	150.0 mm	Natural	10 mL macrotips		3x 100 / bag	N/A	312.10
			Empty rack for 24x 10 mL macrotips		4 / pack	N/A	312.10ER

B = bag, E = empty, F = filter, G = gel, L = long, P = low binding, R = rack, S = sterile, T = refill system

Qualitix® selection chart

Chart indicates optimal combination of Qualitix® and Socorex® pipettes.

10 µL Microtips Tip No.

	3090
	3020
	302
	3090
	3020

20 - 100 µL Microtips Tip No.

	3080
	3080

200 µL Microtips Tip No.

	3070
	3290
	328
	3080
	3070
	3080

300 / 350 µL Microtips Tip No.

	308
	308

1000 µL Microtips Tip No.

	3070
	3190
	3090
	319
	3070

2, 5, 10 mL Macrotips Tip No.

	312
	312
	312

Your distributor:

Acura® manual and Acura® electro																				
815					825 / 826 / 926					835 / 936			855 / 956							
Up to 10 µL	10y / 15 / 20 µL	25 to 100 µL	120 to 200 µL	250 to 1000 µL	0,1 - 2 µL	0,5 - 10 µL	1 - 10 µL (Y)	2 - 20 µL	5 - 50 µL	10 - 100 µL	20 - 200 µL	100 - 1000 µL	0,2 - 2 mL	0,5 - 5 mL	1 - 10 mL	0,5 - 10 µL	5 - 50 µL	10 - 100 µL	20 - 200 µL	40 - 350 µL
X					X	X										X				
X					X	X														
X ²⁾						X ²⁾														
X					X	X										X				
X					X	X														
	X					X	X	X ¹⁾												
		X						X	X	X ¹⁾										
	X	X	X			X	X	X	X	X							X	X	X	
	X	X	X			X	X	X	X	X							X	X	X	
			X					X	X	X	X									
		X	X								X								X	
		X									X									
																				X
																				X
		X										X								
		X										X								
		X										X								
		X										X								
													X							
														X						
															X					

¹⁾ Up to maximal tip volume ²⁾ > 2 µL only



Chemin de Champ-Colomb 7a
 1024 Ecublens / Lausanne
 Switzerland
 Tel. +41 (0)21 651 6000
 socorex@socorex.com
 www.socorex.com



XF.QTX.E - B0917