

# all-glass syringes

## Dosys™ all-glass 155

The reusable borosilicate glass syringes stand for superior chemical and heat shock resistance in various application fields. A precious, low cost alternative to disposable plastic syringes.

*Well suited for  
viscous liquids,  
corrosive reagents,  
dyes, etc.*

### Product advantages

- Precision-machined plunger and barrel
- Excellent fitting and tightness
- Permanent, high visibility graduations
- Sterilisation up to 160°C / 320°F
- Very easy to disassemble and clean



*1 to 50 mL*

### Ordering information - Instruments

Volume	Division	Packaging	Glass Luer nozzle Cat. No.	Metal Luer Lock nozzle, Cat. No.
0.1 - 1 mL	0.05 mL	3 / pk	155.0301	155.0501
0.5 - 2 mL	0.1 mL	3 / pk		155.0502
0.2 - 5 mL	0.2 mL	3 / pk	155.0305	155.0505
1 - 10 mL	0.2 mL	3 / pk	155.0310	155.0510
1 - 20 mL	1 mL	2 / pk	155.0320	155.0520
1 - 50 mL	2 mL	1 / pk		155.0550

### Ordering information - Accessories

#### Dispensing cannulas

Blunt end, superior stainless steel, Luer Lock, autocl.

Description (Ø x l)	Packaging	Cat. No.
0.8 x 40 mm	12 / pk	370.0840
1.2 x 50 mm	1 / pk	370.1250
2.2 x 100 mm	1 / pk	370.22100
2.2 x 150 mm	1 / pk	370.22150

#### Barrel cleaning brushes

Heavy duty nylon, tufted end and plastic handle

Brushes for syringes	Packaging	Cat. No.
1 mL	6 / pk	1.170.301
2, 5 mL	6 / pk	1.170.305
10, 20, 50 mL	3 / pk	1.170.320
Set of two each	6 / pk	1.170.300

For information on our large choice of needles and drenchers or for an overview of the complete Socorex product range, visit [www.socorex.com](http://www.socorex.com)