

Acurex™ compact
more performance – less bulk

low profile dispenser

Superior quality, sturdy dispensers intended for high performance and safe reagent handling. Compact shape, ideal for refrigerator storage. High temperature resistance makes it suitable for water bath and autoclaving.

The Acurex™ line features:

- Integrated dosing mechanism protected inside reservoir
- Simple, robust construction
- Retractable graduated column reduces instrument height
- Fluid path made of inert materials excluding any metal
- Choice of four reservoir sizes, two glass qualities
- Fully autoclavable at 121°C / 250°F

The Acurex™ compact models
501 502



SOCOREX
SWISS

