

all-glass syringes

*Well suited for
viscous liquids,
corrosive reagents,
dyes, etc.*

Dosys™ all-glass 155

The reusable borosilicate glass syringes stand for superior chemical, heat and shock resistance in various application fields.

A precious, low cost alternative to disposable plastic syringes for easy infusing, diluting, aspirating and sampling (not gas tight).

Product advantages

- Precision-machined plunger and barrel
- Excellent fitting and liquid tightness
- Permanent, high visibility graduations
- Glass Luer and metal Luer Lock nozzles
- Sterilisation up to 160°C / 320°F
- Very easy to disassemble and clean



** New *
JUMBO SIZES
up to 250 mL*

Ordering information - Syringes

Volume	Division	Packaging	Glass Luer nozzle Cat. No.	Metal Luer Lock nozzle, Cat. No.
0.1 - 1 mL	0.05 mL	3 / pk	155.0301	155.0501
0.5 - 2 mL	0.1 mL	3 / pk		155.0502
0.2 - 5 mL	0.2 mL	3 / pk	155.0305	155.0505
1 - 10 mL	0.2 mL	3 / pk	155.0310	155.0510
1 - 20 mL	1 mL	2 / pk	155.0320	155.0520
1 - 30 mL	2 mL	2 / pk	155.0330	155.0530
1 - 50 mL	2 mL	1 / pk		155.0550
10 - 100 mL	10 mL	1 / pk	155.03100	155.05100
10 - 150 mL	10 mL	1 / pk		155.05150
10 - 200 mL	10 mL	1 / pk		155.05200
10 - 250 mL	10 mL	1 / pk		155.05250

Dispensing cannulas

Blunt end, superior stainless steel, Luer Lock, autocl.

Description (Ø x l)	Packaging	Cat. No.
0.8 x 40 mm	12 / pk	370.0840
1.2 x 50 mm	1 / pk	370.1250
2.2 x 100 mm	1 / pk	370.22100
2.2 x 150 mm	1 / pk	370.22150

Barrel cleaning brushes

Heavy duty nylon, tufted end and plastic handle

Brushes for syringes	Packaging	Cat. No.
1 mL	6 / pk	1.170.301
2, 5 mL	6 / pk	1.170.305
From 10 to 250 mL	3 / pk	1.170.320
Set of two each	6 / pk	1.170.300

For information on our large choice of needles and drenchers or for an overview of the complete Socorex product range, visit www.socorex.com